

Auxiliary Materials Miscellaneous Products

Name	Viscosity	Thinner	Solids Content approx.	Porcelain	Bone China	Vitreous China	Earthenware	Glass	Notes
Nr. 33	Varnish 85-95 / 6 mm 1)	V 102	10%					●	varnish / liquid crystal adhesive
Nr. 229/1	Fixativ 200-600 mPas 2)	distilled water	2%	●	●	●	●		varnish / for underglaze decals
Nr. 233	Additive	-	-	●	●	●	●	●	softener, additive to the steeping water
Nr. 239	Refreshing Agent 40-45 / 4 mm 1)	-	2%	●	●	●	●	●	re-dissolves the organic binders of a decal and provides the necessary flexibility for correct application

1) Flow time in seconds measured in a DIN cup, temperature: 20°C/68°F

2) Viscosity Thermo Haake RotoVisco 1, Temperature: 20°C/68°F